

## TENDER NAME: OPEN FRAMEWORK AGREEMENT FOR DIABETES

SN	SRM CODE	ITEM NO	ITEM DESCRIPTION	UNIT	QTY	GROUP
1	4111620100600	MOP120005	GLUCOSE MONITOR CONTINUOUS	EACH	1	GROUP 4
2	4200020001100	MOP120005-01	SENSOR RELATED TO GLUCOSE MONITOR CONTINUOUS	EACH	1	GROUP 4
3	4200020001300	MOP120005-02	TRANSMITTER RELATED TO GLUCOSE MONITOR CONTINUOUS	EACH	1	GROUP 4
4	4217000001200	MGE50006	GLUCOMETER	EACH	1	GROUP 1
5	4000000052000	MGE50006-001	GLUCOSE STRIPS RELATED TO GLUCOMETER	EACH	800	GROUP 1
6	4000000052100	MGE50006-002	LANCET RELATED TO GLUCOMETER	EACH	800	GROUP 1
7	4015150500300	MOP120015	PUMP INSULIN MICRO	EACH	1	GROUP 2
8	4200020002100	MOP120015-002	INFUSION SET 6MM RELATED TO PUMP INSULIN MICRO	EACH	360	GROUP 2
9	4200020002200	MOP120015-003	INFUSION SET 9MM RELATED TO PUMP INSULIN MICRO	EACH	360	GROUP 2
10	4200000266700	MOP120015-001	RESERVOIR 2ML-3ML RELATED TO PUMP INSULIN MICRO	EACH	360	GROUP 2
11	4200020002300	MOP120015-004	PATCH RELATED TO PUMP INSULIN MICRO	EACH	360	GROUP 2
12	4200020002400	MOP120015-005	PUMP BASE RELATED TO PUMP INSULIN MICRO	EACH	360	GROUP 2
13	4219261000100	MOP120006	PUMP INSULIN CLOSED LOOP	EACH	1	GROUP 5
14	4200020001400	MOP120006-01	INFUSION SET 6 MM 60-80 CM 90 DGREES INSERTION RELATED TO PUMP INSULIN CLOSED LOOP	EACH	360	GROUP 5
15	4200020001500	MOP120006-02	INFUSION SET 6 MM 110-120 CM 90 DGREES INSERTION RELATED TO PUMP INSULIN CLOSED LOOP	EACH	360	GROUP 5
16	4200020001600	MOP120006-03	INFUSION SET 9 MM 60-80 CM 90 DEGREES INSERTION RELATED TO PUMP INSULIN CLOSED LOOP	EACH	360	GROUP 5
17	4200020001700	MOP120006-04	INFUSION SET 9 MM 110-120 CM 90 DGREES INSERTION RELATED TO PUMP INSULIN CLOSED LOOP	EACH	360	GROUP 5
18	4200020001800	MOP120006-05	INFUSION SET 13 MM 110-120 CM 30-45 DEGREES INSERTION RELATED TO PUMP INSULIN CLOSED LOOP	EACH	360	GROUP 5
19	4200020001900	MOP120006-06	INFUSION SET 13 MM 60-80 CM 30-45 DEGREES INSERTION RELATED TO PUMP INSULIN CLOSED LOOP	EACH	360	GROUP 5
20	4200020002000	MOP120006-07	RESEVOIR 2 ML RELATED TO PUMP INSULIN CLOSED LOOP	EACH	360	GROUP 5
21	4000000049600	MOP120006-08	INSULIN RESERVOIR 3 ML RELATED TO PUMP INSULIN CLOSED LOOP	EACH	360	GROUP 5
22	4000000053500	MOP120006-09	SENSOR INSULIN PUMP RELATED TO PUMP INSULIN CLOSED LOOP	EACH	1	GROUP 5
23	4218190006900	MGE50028	PUMP INSULIN	EACH	1	GROUP 3
24	4200020000400	MGE50028-001	INFUSION SET 6 MM 60-80 CM 90 DEGREES INSERTION RELATED TO PUMP INSULIN	EACH	360	GROUP 3
25	4200020000500	MGE50028-002	INFUSION SET 6 MM 110-120 CM 90 DEGREES INSERTION RELATED TO PUMP INSULIN	EACH	360	GROUP 3
26	4200020000600	MGE50028-003	INFUSION SET 9 MM 60-80 CM 90 DEGREES INSERTION RELATED TO PUMP INSULIN	EACH	360	GROUP 3
27	4200020000700	MGE50028-004	INFUSION SET 9 MM 110-120 CM 90 DEGREES INSERTION RELATED TO PUMP INSULIN	EACH	360	GROUP 3
28	4200020000800	MGE50028-005	INFUSION SET 13 MM 110-120 CM 30 DEGREES INSERTION RELATED TO PUMP INSULIN	EACH	360	GROUP 3
29	4200020000900	MGE50028-006	INFUSION SET 13 MM 60-80 CM 30 DEGREES INSERTION RELATED TO PUMP INSULIN	EACH	360	GROUP 3
30	4200020001000	MGE50028-007	RESEVOIR 2 ML RELATED TO PUMP INSULIN	EACH	360	GROUP 3
31	4000000052700	MGE50028-008	INSULIN RESERVOIR 3 ML OR EQUIVALENT RELATED TO PUMP INSULIN	EACH	360	GROUP 3
32	4111620100600	MOP120005	GLUCOSE MONITOR CONTINUOUS	EACH	1	MOP120005
33	4200020001100	MOP120005-001	SENSOR RELATED TO GLUCOSE MONITOR CONTINUOUS	EACH	1	MOP120005
34	4200020001200	MOP120005-002	INSULIN SMART PEN RELATED TO GLUCOSE MONITOR CONTINUOUS	EACH	1	MOP120005

## TENDER NAME: OPEN FRAMEWORK AGREEMENT FOR DIABETES

SN	SRM CODE	ITEM NO	ITEM DESCRIPTION	UNIT	QTY	GROUP
35	4200020001300	MOP120005-003	TRANSMITTER RELATED TO GLUCOSE MONITOR CONTINUOUS	EACH	1	MOP120005
36	4217000001200	MGE50006	GLUCOMETER	EACH	1	MGE50006
37	4200020000000	MGE50006-003	CARRYING CASE RELATED TO GLUCOMETER	EACH	1	MGE50006
38	4200020000100	MGE50006-004	QUALITY CONTROL TEST KIT LEVEL 1 RELATED TO GLUCOMETER	EACH	1	MGE50006
39	4200020000200	MGE50006-005	QUALITY CONTROL TEST KIT LEVEL 2 RELATED TO GLUCOMETER	EACH	1	MGE50006
40	4200020000300	MGE50006-006	QUALITY CONTROL TEST KIT LEVEL 3 RELATED TO GLUCOMETER	EACH	1	MGE50006
41	4000000052000	MGE50006-001	GLUCOSE STRIPS RELATED TO GLUCOMETER	EACH	1	MGE50006
42	4000000052100	MGE50006-002	LANCET RELATED TO GLUCOMETER	EACH	1	MGE50006
43	4015150500300	MOP120015	PUMP INSULIN MICRO	EACH	1	MOP120015
44	4200000266700	MOP120015-001	RESERVOIR 2ML-3ML RELATED TO PUMP INSULIN MICRO	EACH	1	MOP120015
45	4200020002100	MOP120015-002	INFUSION SET 6MM RELATED TO PUMP INSULIN MICRO	EACH	1	MOP120015
46	4200020002200	MOP120015-003	INFUSION SET 9MM RELATED TO PUMP INSULIN MICRO	EACH	1	MOP120015
47	4200020002300	MOP120015-004	PATCH RELATED TO PUMP INSULIN MICRO	EACH	1	MOP120015
48	4200020002400	MOP120015-005	PUMP BASE RELATED TO PUMP INSULIN MICRO	EACH	1	MOP120015
49	4219261000100	MOP120006	PUMP INSULIN CLOSED LOOP	EACH	1	MOP120006
50	4200020001400	MOP120006-001	INFUSION SET 6 MM 60-80 CM 90 DGREES INSERTION RELATED TO PUMP INSULIN CLOSED LOOP	EACH	1	MOP120006
51	4200020001500	MOP120006-002	INFUSION SET 6 MM 110-120 CM 90 DGREES INSERTION RELATED TO PUMP INSULIN CLOSED LOOP	EACH	1	MOP120006
52	4200020001600	MOP120006-003	INFUSION SET 9 MM 60-80 CM 90 DEGREES INSERTION RELATED TO PUMP INSULIN CLOSED LOOP	EACH	1	MOP120006
53	4200020001700	MOP120006-004	INFUSION SET 9 MM 110-120 CM 90 DGREES INSERTION RELATED TO PUMP INSULIN CLOSED LOOP	EACH	1	MOP120006
54	4200020001800	MOP120006-005	INFUSION SET 13 MM 110-120 CM 30-45 DEGREES INSERTION RELATED TO PUMP INSULIN CLOSED LOOP	EACH	1	MOP120006
55	4200020001900	MOP120006-006	INFUSION SET 13 MM 60-80 CM 30-45 DEGREES INSERTION RELATED TO PUMP INSULIN CLOSED LOOP	EACH	1	MOP120006
56	4200020002000	MOP120006-007	RESEVOIR 2 ML RELATED TO PUMP INSULIN CLOSED LOOP	EACH	1	MOP120006
57	4000000049600	MOP120006-008	INSULIN RESERVOIR 3 ML RELATED TO PUMP INSULIN CLOSED LOOP	EACH	1	MOP120006
58	4000000053500	MOP120006-009	SENSOR INSULIN PUMP RELATED TO PUMP INSULIN CLOSED LOOP	EACH	1	MOP120006
59	4218190006900	MGE50028	PUMP INSULIN	EACH	1	MGE50028
60	4200020000400	MGE50028-001	INFUSION SET 6 MM 60-80 CM 90 DEGREES INSERTION RELATED TO PUMP INSULIN	EACH	1	MGE50028
61	4200020000500	MGE50028-002	INFUSION SET 6 MM 110-120 CM 90 DEGREES INSERTION RELATED TO PUMP INSULIN	EACH	1	MGE50028
62	4200020000600	MGE50028-003	INFUSION SET 9 MM 60-80 CM 90 DEGREES INSERTION RELATED TO PUMP INSULIN	EACH	1	MGE50028
63	4200020000700	MGE50028-004	INFUSION SET 9 MM 110-120 CM 90 DEGREES INSERTION RELATED TO PUMP INSULIN	EACH	1	MGE50028
64	4200020000800	MGE50028-005	INFUSION SET 13 MM 110-120 CM 30 DEGREES INSERTION RELATED TO PUMP INSULIN	EACH	1	MGE50028
65	4200020000900	MGE50028-006	INFUSION SET 13 MM 60-80 CM 30 DEGREES INSERTION RELATED TO PUMP INSULIN	EACH	1	MGE50028
66	4200020001000	MGE50028-007	RESEVOIR 2 ML RELATED TO PUMP INSULIN	EACH	1	MGE50028
67	4000000052700	MGE50028-008	INSULIN RESERVOIR 3 ML OR EQUIVALENT RELATED TO PUMP INSULIN	EACH	1	MGE50028